

# **Model: Collection tray of sheet steel galvanised for 2 x 200 l-drums Order number: 1.39.202**

Correct use: Storage equipment and accessories



dimensions: approx. 1200 x 800 x 360 mm (L x W x H) collection volume: 215 l material strength 3 mm, with galvanised removable gratings with supporting feet, 100 mm underclearance suitable for flammable liquids F+, F and water polluting liquids WGK 1 - 3 Testings: with declaration of conformity (ÜHP) according to StawaR Tested by TÜV

Category: Collecting trays Outer height [mm]: 360, Outer width [mm]: 1.200, Outer depth [mm]: 800, Gross weight [kg]: 72,00, Net weight [kg]: 72,00,