

**Model: Collection tray of sheet steel galvanised
for 2 x 200 l-drums
Order number: 1.39.204**

Correct use: Storage equipment and accessories



- square - dimensions: approx. 1200 x 1200 x 285 mm (LxWxH) collection volume: 216 l
material strength 3 mm, with galvanised removable gratings with supporting feet, 100 mm
underclearance suitable for flammable liquids F+, F and water polluting liquids WGK 1 - 3
Testings: with declaration of conformity (ÜHP) according to StawaR Tested by TÜV

Category:	Collecting trays
Outer width [mm]:	1.200
Outer depth [mm]:	1.200
Outer height [mm]:	285
Net weight [kg]:	91,00