

## **Model: Collection tray of sheet steel galvanised for 2 x 200 l-drums**

**Order number: 1.39.204**

Correct use: Storage equipment and accessories



- square - dimensions: approx. 1200 x 1200 x 285 mm (LxWxH) collection volume: 216 l  
material strength 3 mm, with galvanised removable gratings with supporting feet, 100 mm  
underclearance suitable for flammable liquids F+, F and water polluting liquids WGK 1 - 3  
Testings: with declaration of conformity (ÜHP) according to StawaR Tested by TÜV

Category: Collecting trays Outer height [mm]: 285, Outer width [mm]: 1.200, Outer depth [mm]:  
1.200, Gross weight [kg]: 91,00, Net weight [kg]: 91,00,