


**Model: Collection tray of sheet steel painted for
10 x 200 l-drums**
Order number: 1.40.501

Correct use: Storage equipment and accessories



dimensions: approx. 3280 x 1300 x 200 (L x W x H) in mm, collection volume: 280 l colour RAL 5012 blue, material strength 3 mm, with galvanised, removable grating (3 pieces), fork stirrups 100 mm underclearance, suitable for flammable liquids F+, F and water polluting liquids WGK 1-3
Testings: with declaration of conformity (ÜHP) according to StawaR Tested by TÜV

Category:	Collecting trays
Outer width [mm]:	3.280
Outer depth [mm]:	1.300
Outer height [mm]:	200
Net weight [kg]:	242,00