

Model: Collection tray of sheet steel painted for 20 x 200 l-drums

Order number: 1.40.502

Correct use: Storage equipment and accessories



dimensions: approx. 3880 x 1300 x 200, (L x W x H) in mm, collection volume: 335 l colour RAL 5012 blue, material strength 3 mm, with galvanised, removable grating (3 pieces), fork stirrups 100 mm underclearance, suitable for flammable liquids F+ F and water polluting liquids WGK 1-3
Testings: with declaration of conformity (ÜHP) according to StawaR Tested by TÜV

Category: Collecting trays Outer height [mm]: 200, Outer width [mm]: 3.880, Outer depth [mm]: 1.300, Gross weight [kg]: 277,00, Net weight [kg]: 277,00,