

# **Model: Collection tray of sheet steel galvanised for 10 x 200 l-drum Order number: 1.40.601**

Correct use: Storage equipment and accessories



dimensions: approx. 3280 x 1300 x 200, (L x W x H) in mm, collection volume: 280 l material strength 3 mm, with galvanised, removable grating (3 pieces), fork stirrups 100 mm underclearance, suitable for flammable liquids F+ F and water polluting liquids WGK 1-3 Testings: with declaration of conformity (ÜHP) according to StawaR Tested by TÜV

Category: Collecting trays Outer height [mm]: 200, Outer width [mm]: 3.280, Outer depth [mm]: 1.300, Gross weight [kg]: 262,00, Net weight [kg]: 260,00,