

**Model:**  
**Order number: 1.41.202-1**

Correct use: Storage equipment and accessories



Approx. dim.: 1200 x 800 x 360 (WxDxH in mm) Collection volume: 200 l Colour RAL 5010 blue, material strength 3 mm, with galvanised removable grating, skids with 100 mm underclearance, Load: 500 kg Suitable for flammable liquids F+, F und for water polluting liquids WGK 1-3 Testings: Declaration of conformity (ÜHP) according to StawaR,

Category:	Collecting trays
Outer width [mm]:	1.200
Outer depth [mm]:	800
Outer height [mm]:	300
Net weight [kg]:	70,00