

Model:

Order number: 1.41.202-1

Correct use: Storage equipment and accessories



Approx. dim.: 1200 x 800 x 360 (WxDxH in mm) Collection volume: 200 l Colour RAL 5010 blue, material strength 3 mm, with galvanised removable grating, skids with 100 mm underclearance, Load: 500 kg Suitable for flammable liquids F+, F und for water polluting liquids WGK 1-3 Testings: Declaration of conformity (ÜHP) according to StawaR,

Category: Collecting trays Outer height [mm]: 300, Outer width [mm]: 1.200, Outer depth [mm]: 800, Gross weight [kg]: 70,00, Net weight [kg]: 70,00,