

Model: Collection tray of sheet steel painted for 4 x 200 l-drum

Order number: 1.41.204

Correct use: Storage equipment and accessories



Dimensions: appr. 2400 x 800 x 250 (WxDxH in mm) Collection volume: 210 l Colour RAL 5012 blue, material strength 3 mm, with galvanised removable grating, skids with 100 mm underclearance, Load: 1000 kg Suitable for flammable liquids F+, F und for water polluting liquids WGK 1-3 Testings: Declaration of conformity (ÜHP) according to StawaR, TÜV-tested

Category: Collecting trays Outer height [mm]: 250, Outer width [mm]: 2.400, Outer depth [mm]: 800, Gross weight [kg]: 113,00, Net weight [kg]: 113,00,