

Model: Collection tray of sheet steel galvanised for 1 x 200 l-drum

Order number: 1.41.301

Correct use: Storage equipment and accessories



Approx. dim.: 800 x 800 x 465 (WxDxH in mm) Collection volume: 210 l Material strength 3 mm, with galvanised removable grating, skids with 100 mm underclearance, Load: 250 kg Suitable for flammable liquids F+, F and for water polluting liquids WGK 1-3 Testings: Declaration of conformity (ÜHP) according to StawaR, TÜV-tested

Category:	Collecting trays
Outer width [mm]:	800
Outer depth [mm]:	800
Outer height [mm]:	465
Net weight [kg]:	63,00