

Model: Collection tray of sheet steel galvanised for 1 x 200 l-drum Order number: 1.41.301

Correct use: Storage equipment and accessories



Approx. dim.: 800 x 800 x 465 (WxDxH in mm) Collection volume: 210 l Material strength 3 mm, with galvanised removable grating, skids with 100 mm underclearance, Load: 250 kg Suitable for flammable liquids F+, F and for water polluting liquids WGK 1-3 Testings: Declaration of conformity (ÜHP) according to StawaR, TÜV-tested

Category: Collecting trays Outer height [mm]: 465, Outer width [mm]: 800, Outer depth [mm]: 800, Gross weight [kg]: 63,00, Net weight [kg]: 63,00,