

### **Model:**

**Order number: 1.41.302-1**

Correct use: Storage equipment and accessories



Approx. dim.: 1200 x 800 x 360 (WxDxH in mm) Collection volume: 210 l Material strength 3 mm, with galvanised removable grating, skids with 100 mm underclearance, Load: 500 kg Suitable for flammable liquids F+, F und for water polluting liquids WGK 1-3 Testings: Declaration of conformity (ÜHP) according to StawaR

Category: Collecting trays Outer height [mm]: 360, Outer width [mm]: 1.200, Outer depth [mm]: 800, Gross weight [kg]: 75,00, Net weight [kg]: 75,00,