

Model: Collection tray of sheet steel galvanised for 4 x 200 l-drum

Order number: 1.41.304

Correct use: Storage equipment and accessories



Approx. dim.: 2400 x 800 x 250 (WxDxH in mm) Collection volume: 210 l Material strength 3 mm, with galvanised removable grating, skids with 100 mm underclearance, Load: 1000 kg Suitable for flammable liquids F+, F und for water polluting liquids WGK 1-3 Testings: Declaration of conformity (ÜHP) according to StawaR, TÜV-tested

Category: Collecting trays Outer height [mm]: 250, Outer width [mm]: 2.400, Outer depth [mm]: 800, Gross weight [kg]: 120,00, Net weight [kg]: 120,00,