

**Model: Hazardous substance shelving extension
for art. 1.46.530*
Order number: 1.47.530**

Correct use: Storage equipment and accessories



with 4 x Basin shelf Dim.: 1000 x 600 x 2000 mm (W x D x H) Collection volume: 4 x 30 litre
Suitable for flammable liquids F+, F and water polluting liquids WGK 1-3 Galvanised shelving
uprights and basin shelves, height adjustable with 25 mm . Load bearing capacity per shelf 150 kg
for uniformly distributed load Shelf stick system to self installation testing Collection trays with
declaration of conformity (ÜHP) according to StawaR (Steel collection tray regulations), DLG-fokus-
test, TÜV-tested

Category:	Shelf systems for hazardous substances
Outer width [mm]:	1.000
Outer depth [mm]:	600
Outer height [mm]:	2.000
Net weight [kg]:	114,00