

# **Model: Hazardous substance shelving extension for art. 1.46.530\***

## **Order number: 1.47.530**

Correct use: Storage equipment and accessories



with 4 x Basin shelf Dim.: 1000 x 600 x 2000 mm (W x D x H) Collection volume: 4 x 30 litre  
Suitable for flammable liquids F+, F and water polluting liquids WGK 1-3 Galvanised shelving  
uprights and basin shelves, height adjustable with 25 mm . Load bearing capacity per shelf 150 kg  
for uniformly distributed load Shelf stick system to self installation testing Collection trays with  
declaration of conformity (ÜHP) according to StawaR (Steel collection tray regulations), DLG-fokus-  
test, TÜV-tested

Category: Shelf systems for hazardous substances Outer height [mm]: 2.000, Outer width [mm]:  
1.000, Outer depth [mm]: 600, Gross weight [kg]: 114,00, Net weight [kg]: 114,00,