

Model: Hazardous substance depot, galvanised, with two wing doors

Order number: 1.49.212

Correct use: Storage equipment and accessories



galvanised with 2 wings doors for 2 x 200 l-drums Dims.: 1437 x 1015 x 1614 mm (L x W x H)
Collection vol.: 210 l painted external surface blue RAL 5012 with natural air exchange collection
tray with galvanised grating Roof with pneumatic springs, lockable fork stirrups with 100 mm
underclearance for flammable liquids F+, F and water polluting liquids in water pollution class WGK
1-3. Testing: With general national technical approval Z-38.5-101 (DIBt Berlin), TÜV-tested

Category: Hazardous substances depot Outer height [mm]: 1.614, Outer width [mm]: 1.437, Outer
depth [mm]: 1.015, Gross weight [kg]: 236,00, Net weight [kg]: 236,00,