

**Model: Hazardous substance depot, galvanised,
with two wing doors**
Order number: 1.49.214

Correct use: Storage equipment and accessories



galvanised with 2 wings doors for 4 x 200 l-drums, Dims.: 1437 x 1500 x 1557 mm (L x W x H)
Collection vol.: 220 l with natural air exchange collection tray with galvanised grating Roof with
pneumatic springs, lockable fork stirrups with 100 mm underclearance for flammable liquids F+, F
and water polluting liquids in water pollution class WGK 1-3. Testing: With general national
technical approval Z-38.5-101 (DIBt Berlin), TÜV-tested

Category:	Hazardous substances depot
Outer width [mm]:	1.437
Outer depth [mm]:	1.500
Outer height [mm]:	1.557
Net weight [kg]:	302,00