

Model: Hazardous substance depot, galvanised, with two wing doors

Order number: 1.49.214

Correct use: Storage equipment and accessories



galvanised with 2 wings doors for 4 x 200 l-drums, Dims.: 1437 x 1500 x 1557 mm (L x W x H)
Collection vol.: 220 l with natural air exchange collection tray with galvanised grating Roof with
pneumatic springs, lockable fork stirrups with 100 mm underclearance for flammable liquids F+, F
and water polluting liquids in water pollution class WGK 1-3. Testing: With general national
technical approval Z-38.5-101 (DIBt Berlin), TÜV-tested

Category: Hazardous substances depot Outer height [mm]: 1.557, Outer width [mm]: 1.437, Outer
depth [mm]: 1.500, Gross weight [kg]: 302,00, Net weight [kg]: 302,00,