

**Model: Collection tray, PE, for 2 x 200 l -drum,  
upright**  
**Order number: 1.61.002**

Correct use: Storage equipment and accessories



Dim.: appr. 1270 x 870 x 400 (W x D x H in mm), collection volume: 225 l load bearing capacity: 900 kg low tare weight: 62 kg plastic structure made from extruded laminate (PE-HD), with stable PE perforated plate as standing level accessible due to 100 mm base design. suitable for water polluting, non-flammable liquids WGK 1-3, Tested: General national technical approval of DIBT-Berlin, MPA-tested, TÜV-tested

Category:	Collecting trays
Outer width [mm]:	1.270
Outer depth [mm]:	870
Outer height [mm]:	400
Net weight [kg]:	62,00