

Model: Collection tray, PE, for 2 x 200 l -drum, upright

Order number: 1.61.002

Correct use: Storage equipment and accessories



Dim.: appr. 1270 x 870 x 400 (W x D x H in mm), collection volume: 225 l load bearing capacity: 900 kg low tare weight: 62 kg plastic structure made from extruded laminate (PE-HD), with stable PE perforated plate as standing level accessible due to 100 mm base design. suitable for water polluting, non-flammable liquids WGK 1-3, Tested: General national technical approval of DIBT-Berlin, MPA-tested, TÜV-tested

Category: Collecting trays Outer height [mm]: 400, Outer width [mm]: 1.270, Outer depth [mm]: 870, Gross weight [kg]: 64,00, Net weight [kg]: 62,00,