

**Model: Collection tray, PE, for 3 x 200 l-drum,
upright**
Order number: 1.61.003

Correct use: Storage equipment and accessories



Dim.: appr. 1870 x 870 x 310 (W x D x H in mm), collection volume: 225 l load bearing capacity: 1350 kg low tare weight: 77 kg plastic structure made from extruded laminate (PE-HD), with stable PE perforated plate as standing level accessible due to 100 mm base design. suitable for water polluting, non-flammable liquids WGK 1-3, Tested: General national technical approval of DIBT-Berlin, MPA-tested, TÜV-tested

Category:	Collecting trays
Outer width [mm]:	1.870
Outer depth [mm]:	870
Outer height [mm]:	310
Net weight [kg]:	77,00