

Model: Collection tray, PE, for 3 x 200 l-drum, upright

Order number: 1.61.003

Correct use: Storage equipment and accessories



Dim.: appr. 1870 x 870 x 310 (W x D x H in mm), collection volume: 225 l load bearing capacity: 1350 kg low tare weight: 77 kg plastic structure made from extruded laminate (PE-HD), with stable PE perforated plate as standing level accessible due to 100 mm base design. suitable for water polluting, non-flammable liquids WGK 1-3, Tested: General national technical approval of DIBT-Berlin, MPA-tested, TÜV-tested

Category: Collecting trays Outer height [mm]: 310, Outer width [mm]: 1.870, Outer depth [mm]: 870, Gross weight [kg]: 79,00, Net weight [kg]: 77,00,