

# **Model: Collection tray, PE, for 4 x 200 l-drum, upright**

## **Order number: 1.61.004**

Correct use: Storage equipment and accessories



Dim.: appr. 2470 x 870 x 270 (W x D x H in mm), collection volume: 225 l load bearing capacity: 1800 kg low tare weight: 89 kg plastic structure made from extruded laminate (PE-HD), with stable PE perforated plate as standing level accessible due to 100 mm base design. suitable for water polluting, non-flammable liquids WGK 1-3, Tested: General national technical approval of DIBT-Berlin, MPA-tested, TÜV-tested

Category: Collecting trays Outer height [mm]: 270, Outer width [mm]: 2.470, Outer depth [mm]: 870, Gross weight [kg]: 90,00, Net weight [kg]: 89,00,