

**Model: Collection tray, PE, with grating, collection  
vol: 40 l  
Order number: 1.73.040**

Correct use: Storage equipment and accessories



for small containers dim.: approx. 1300 x 620 x 80 (L x W x H in mm), load bearing capacity : 50 kg, low tare weight: 13 kg, with galvanised steel grating suitable for water polluting liquids WGK 1 - 3. testing: with general national technical approval Z-40.22-388 of the DIBt-Berlin

Category:	Collecting trays
Outer width [mm]:	1.300
Outer depth [mm]:	620
Outer height [mm]:	80
Net weight [kg]:	13,00