

Model: Collection tray, PE, with grating, collection vol: 40 l Order number: 1.73.040

Correct use: Storage equipment and accessories



for small containers dim.: approx. 1300 x 620 x 80 (L x W x H in mm), load bearing capacity : 50 kg, low tare weight: 13 kg, with galvanised steel grating suitable for water polluting liquids WGK 1 - 3. testing: with general national technical approval Z-40.22-388 of the DIBt-Berlin

Category: Collecting trays Outer height [mm]: 80, Outer width [mm]: 1.300, Outer depth [mm]: 620, Gross weight [kg]: 13,00, Net weight [kg]: 13,00,