

Model: Chemical absorber Extreme tear strenght, 2 rolls

Order number: 11.51.300

Correct use: Absorber and binders



non-woven fabric binders made from polypropylene micro fibres in plastic bag dimensions: approx. 0.40 x 44 m can be individually cut to required size maximum absorption capacity per VE: 163 l
Testing: tested by MPA NRW, Dortmund and Hygiene-Institut, Gelsenkirchen

Category: Nonwovens fabrics made of micro-fibres Gross weight [kg]: 10,50, Net weight [kg]: 8,60,