

# **Model: Emergency kit in transport box**

## **Absorbency 75 l**

### **Order number: 11.55.000**

Correct use: Absorber and binders



lockable transport box with rolls and handle, made from highly stable PVC Dim.: approx. 80 x 38 x 43 cm (L x W x H) with universal absorbers - extra absorbent - for industrial fluids: 35 absorber cloths, 3 absorber socks (diam.x L ca. 7 cm x 1,20 m) 2 pads (non-woven binders made from polypropylene micro fibres) plus 10 cloths, 1 pair of gloves, 3 bags Testing: absorbent tested by MPA NRW, Dortmund and Hygiene-Institut, Gelsenkirchen

Category: Absorber emergency kits Gross weight [kg]: 13,00, Net weight [kg]: 9,70,