

**Model: Additional pull-out tray BENCH S**  
**Order number: 29-20661-119**

Correct use: Passive storage



- With the additional pull-out tray mounted, the volume of the largest individual container changes to 2,5l and the max. load capacity to 35 kg.
- Setting heights: bottom: 270 mm, top: 270 mm,
- powder-coated sheet steel RAL 7035 light grey

Category:	Equipment
Outer width [mm]:	425
Outer depth [mm]:	433
Outer height [mm]:	350
Net weight [kg]:	11,00