

Model: Additional pull-out tray BENCH S

Order number: 29-20661-119

Correct use: Passive storage



- With the additional pull-out tray mounted, the volume of the largest individual container changes to 2,5l and the max. load capacity to 35 kg.
- Setting heights: bottom: 270 mm, top: 270 mm,
- powder-coated sheet steel RAL 7035 light grey

Category: Equipment Outer height [mm]: 350, Outer width [mm]: 425, Outer depth [mm]: 433,
Gross weight [kg]: 11,00, Net weight [kg]: 11,00,