

Model: Locker (door hung left) separately lockable

Order number: 29-50000-036L

Correct use: Passive storage



approx. dimensions: 370 x 427 x 359 (WxDxH) mm powder-coated sheet steel RAL 7035 light-grey with perforated side walls, for optimized ventilation

| | |
|--------------------|-----------|
| Category: | Equipment |
| Outer width [mm]: | 370 |
| Outer depth [mm]: | 427 |
| Outer height [mm]: | 359 |
| Net weight [kg]: | 16,00 |