

Model: Sheet steel tilting station capacity 19 l

Order number: 3.108.88

Correct use: Safety containers and cleaning containers



Consisting of: Tilt can 19 litres and framework rack can made from sheet steel, interior tin plated
Exterior red powder coated Dosing tap from brass Self-closing openings Flame arrester in in filling
opening Automatic pressure / vacuum compansation Tipping frame : 365 x 370 x 600 (L x W x H)
mm Marking prepared according to GefStoffV and BGV A8 Testing: FM-approved

Category:	for liquid hazardous substances
Outer width [mm]:	365
Outer depth [mm]:	370
Outer height [mm]:	600
Net weight [kg]:	6,40