

# **Model: Sheet steel tilting station capacity 19 l**

## **Order number: 3.108.88**

Correct use: Safety containers and cleaning containers



Consisting of: Tilt can 19 litres and framework rack can made from sheet steel, interior tin plated  
Exterior red powder coated Dosing tap from brass Self-closing openings Flame arrester in in filling  
opening Automatic pressure / vacuum compansation Tipping frame : 365 x 370 x 600 (L x W x H)  
mm Marking prepared according to GefStoffV and BGV A8 Testing: FM-approved

Category: for liquid hazardous substances Outer height [mm]: 600, Outer width [mm]: 365, Outer  
depth [mm]: 370, Gross weight [kg]: 6,50, Net weight [kg]: 6,40,