

# **Model: Steel moistener capacity 0,5 l**

## **Order number: 4.100.08**

Correct use: Safety containers and cleaning containers



Dimensions: 120 x 135 ( dia x H in mm) Soaking tray: 75 mm dia, interior tin plated, exterior red powder coated, soaking tray with flame arrester, brass pump tube and PPS, Markings prepared according to GefStoffV and BGV A8. testing: FM-tested

Category: Cleaning containers Outer height [mm]: 135, Gross weight [kg]: 1,00, Net weight [kg]: 0,50,