

Model: Sheet steel moistener capacity 2 l

Order number: 4.102.08

Correct use: Safety containers and cleaning containers



Dimensions: 180 x 190 (dia x H in mm). Soaking tray: 135 mm dia, interior tin plated, exterior red powder coated, soaking tray with flame arrester, brass pump tube and PPS. Marking prepared according to GefStoffV and BGV A8. testing: FM-tested

Category: Cleaning containers Outer height [mm]: 190, Gross weight [kg]: 1,50, Net weight [kg]: 1,00,