

Model: Sheet steel moistener capacity 4 l

Order number: 4.103.08

Correct use: Safety containers and cleaning containers



Dimensions: 180 x 270 (dia x H in mm). Soaking tray: 135 mm dia, interior tin plated, exterior red powder coated, soaking tray with flame arrester, brass pump tube and PPS. Marking prepared according to GefStoffV and BGV A8. testing: FM-tested

Category: Cleaning containers Outer height [mm]: 270, Gross weight [kg]: 2,00, Net weight [kg]: 1,40,