


**Model: Parts basket for dipping tank model
4.271.10 and 4.273.11
Order number: 4.270.03**

Correct use: Safety containers and cleaning containers



Galvanised parts basket for dipping tank model 4.271.10 and 4.273.11 Dim. (WxDxH): approx. 406 x 330 x 102 mm. Offers effective cleaning of multiple parts for every cleaning section.

Category:	Cleaning tanks
Outer width [mm]:	406
Outer depth [mm]:	330
Outer height [mm]:	102
Net weight [kg]:	1,80