


# **Model: Parts basket for dipping tank model 4.271.10 and 4.273.11 Order number: 4.270.03**

Correct use: Safety containers and cleaning containers



Galvanised parts basket for dipping tank model 4.271.10 and 4.273.11 Dim. (WxDxH): approx. 406 x 330 x 102 mm. Offers effective cleaning of multiple parts for every cleaning section.

Category: Cleaning tanks Outer height [mm]: 102, Outer width [mm]: 406, Outer depth [mm]: 330, Gross weight [kg]: 4,00, Net weight [kg]: 1,80,