

**Model: Strainer for dipping tank model 4.271.10
and 4.273.11**
Order number: 4.270.04

Correct use: Safety containers and cleaning containers



Galvanised strainer for dipping tank model 4.271.10 and 4.273.11 Dim. (WxDxH): approx. 584 x 349 x 19 mm. Gathers coarse dirt caused by sedimentation. Offers easy cleaning during solvent exchanging.

Category:	Cleaning tanks
Outer width [mm]:	584
Outer depth [mm]:	349
Outer height [mm]:	19
Net weight [kg]:	0,90