


Model: Strainer for dipping tank model 4.271.10 and 4.273.11

Order number: 4.270.04

Correct use: Safety containers and cleaning containers



Galvanised strainer for dipping tank model 4.271.10 and 4.273.11 Dim. (WxDxH): approx. 584 x 349 x 19 mm. Gathers coarse dirt caused by sedimentation. Offers easy cleaning during solvent exchanging.

Category: Cleaning tanks Outer height [mm]: 19, Outer width [mm]: 584, Outer depth [mm]: 349, Gross weight [kg]: 2,00, Net weight [kg]: 0,90,