


Model: Parts basket for dipping tank model 4.272.20 and 4.273.22 Order number: 4.270.05

Correct use: Safety containers and cleaning containers



Galvanised parts basket for dipping tank model 4.271.20 and 4.273.22 Dim. (WxDxH): approx. 711 x 330 x 102 mm. Offers effective cleaning of multiple parts every cleaning section.

Category: Cleaning tanks Outer height [mm]: 102, Outer width [mm]: 711, Outer depth [mm]: 330, Gross weight [kg]: 5,00, Net weight [kg]: 2,30,