

Model: Pedal operated dipping tank, sheet steel, 48 l
Order number: 4.271.10

Correct use: Safety containers and cleaning containers



capacity: approx. 48 litres outside dimensions: 610 x 496 x 845 (WxDxH) mm basin (inside): 603 x 406 x 222 (WxDxH) mm Basin, lid and frame made from galvanised sheet steel, red coated. "3/4"-drain plug in basin bottom for selective emptying. Slow, independent closing of the lid via pneumatic cylinder. Testing: FM-approved

Category:	Cleaning tanks
Outer width [mm]:	610
Outer depth [mm]:	406
Outer height [mm]:	222
Net weight [kg]:	27,20