

Model: Pedal operated dipping tank, sheet steel, 48 l

Order number: 4.271.10

Correct use: Safety containers and cleaning containers



capacity: approx. 48 litres outside dimensions: 610 x 496 x 845 (WxDxH) mm basin (inside): 603 x 406 x 222 (WxDxH) mm Basin, lid and frame made from galvanised sheet steel, red coated. "3/4"-drain plug in basin bottom for selective emptying. Slow, independent closing of the lid via pneumatic cylinder. Testing: FM-approved

Category: Cleaning tanks Outer height [mm]: 222, Outer width [mm]: 610, Outer depth [mm]: 406, Gross weight [kg]: 30,00, Net weight [kg]: 27,20,