

**Model: Dipping tank, sheet steel 2 x 15 l
chambers**
Order number: 4.271.40

Correct use: Safety containers and cleaning containers



capacity: approx. 2 x 15 litres outside dimensions: 610 x 496 x 845 mm (W x D x H) basin (inside): 2x 298 x 406 x 140 mm(W x D x H) Basin, lid and frame made from galvanised sheet steel, red coated. For cleaning of extremely contaminated components. 1. chamber: pre-cleaning, 2. chamber: final cleaning "3/4""- drain plug in basin bottom for" selective emptying. Testing: FM-approved

Category:	Cleaning tanks
Outer width [mm]:	610
Outer depth [mm]:	406
Outer height [mm]:	222
Net weight [kg]:	28,10