

Model: Dipping tank, sheet steel 2 x 15 l chambers

Order number: 4.271.40

Correct use: Safety containers and cleaning containers



capacity: approx. 2 x 15 litres outside dimensions: 610 x 496 x 845 mm (W x D x H) basin (inside): 2x 298 x 406 x 140 mm(W x D x H) Basin, lid and frame made from galvanised sheet steel, red coated. For cleaning of extremely contaminated components. 1. chamber: pre-cleaning, 2. chamber: final cleaning "3/4"- drain plug in basin bottom for" selective emptying. Testing: FM-approved

Category: Cleaning tanks Outer height [mm]: 222, Outer width [mm]: 610, Outer depth [mm]: 406, Gross weight [kg]: 30,50, Net weight [kg]: 28,10,