

Model: Pedal operated dipping tank, sheet steel, 97 litres

Order number: 4.272.20

Correct use: Safety containers and cleaning containers



capacity: 97 litres outside dimensions: 895 x 496 x 845 (WxDxH) mm basin (inside): 899 x 406 x 305 (WxDxH) mm Basin, lid and frame made from galvanised sheet steel, red coated. "3/4"- drain plug in basin bottom for" selective emptying. Slow, independent closing of lid via pneumatic cylinder. testing:

Category:	Cleaning tanks
Outer width [mm]:	895
Outer depth [mm]:	406
Outer height [mm]:	305
Net weight [kg]:	35,40