

Model: Pedal operated dipping tank, sheet steel, 97 litres Order number: 4.272.20

Correct use: Safety containers and cleaning containers



capacity: 97 litres outside dimensions: 895 x 496 x 845 (WxDxH) mm basin (inside): 899 x 406 x 305 (WxDxH) mm Basin, lid and frame made from galvanised sheet steel, red coated. "3/4"- drain plug in basin bottom for" selective emptying. Slow, independent closing of lid via pneumatic cylinder. testing:

Category: Cleaning tanks Outer height [mm]: 305, Outer width [mm]: 895, Outer depth [mm]: 406, Gross weight [kg]: 37,50, Net weight [kg]: 35,40,