

Model: Dipping tank, sheet steel with chambers 2 x 34 l

Order number: 4.272.60

Correct use: Safety containers and cleaning containers



capacity: approx. 2 x 34 litres Outside dimensions: 895 x 496 x 845 mm (WxDxH) basin (inside): 2x438 x 406 x 222 mm (WxDxH) Basin, lid and frame made from galvanised sheet steel, red coated. For cleaning of extremely contaminated components. 1. chamber: pre-cleaning, 2. chamber: final purification. "3/4"- drain plug in bottom of each chamber" for selective emptying. Testing: FM-approved

Category: Cleaning tanks Outer height [mm]: 305, Outer width [mm]: 895, Outer depth [mm]: 406, Gross weight [kg]: 38,50, Net weight [kg]: 36,30,