

Model: Dipping tank for workbenches sheet steel, 48 l
Order number: 4.273.11

Correct use: Safety containers and cleaning containers



capacity: approx. 48 litres outside dimensions: 610 x 496 x 273 mm (W x D x H) basin (inside): 603 x 406 x 273 mm (W x D x H) For positioning on a workbench. Basin and lid made from galvanised sheet steel, red coated. self-closing lid in case of fire "3/4""- drain plug in basin bottom for selective emptying. Testing: FM-approved

Category:	Cleaning tanks
Outer width [mm]:	610
Outer depth [mm]:	406
Outer height [mm]:	222
Net weight [kg]:	13,30