

# **Model: Dipping tank for workbenches sheet steel, 48 l**

## **Order number: 4.273.11**

Correct use: Safety containers and cleaning containers



capacity: approx. 48 litres outside dimensions: 610 x 496 x 273 mm (W x D x H) basin (inside): 603 x 406 x 273 mm (W x D x H) For positioning on a workbench. Basin and lid made from galvanised sheet steel, red coated. self-closing lid in case of fire "3/4"- drain plug in basin bottom for selective emptying. Testing: FM-approved

Category: Cleaning tanks Outer height [mm]: 222, Outer width [mm]: 610, Outer depth [mm]: 406, Gross weight [kg]: 15,50, Net weight [kg]: 13,30,