

Model: Dipping tank for work benches sheet steel, 97 l

Order number: 4.273.22

Correct use: Safety containers and cleaning containers



capacity: approx. 97 litres outside dimensions: 895 x 496 x 336 mm (W x D x H) basin (inside): 889 x 406 x 336 mm (W x D x H) For positioning on a workbench. Basin, lid and frame made from galvanised sheet steel, red coated. self-closing lid in case of fire "3/4"- drain plug in basin bottom for selective emptying Testing: FM-approved

Category:	Cleaning tanks
Outer width [mm]:	895
Outer depth [mm]:	406
Outer height [mm]:	305
Net weight [kg]:	22,00